

Table of Content

1.	Release note 001	2
1.1	New feature	2
1.2	Feature Adjustments.....	2
1.3	Bug Fixes.....	2
2.	Release note 002	3
2.1	New feature	3
2.2	Feature Adjustments.....	3
2.3	Bug Fixes.....	3
	File records.....	4

1. Release note 001

Date: 2024.06.04

1.1 New feature

- Task creation, Type selection, Task start, Termination, Result viewing
- View historical, Pagination, Deletion
- System setting, Language, Parameter settings
- Reaction condition parameters as shown in the following diagram

Test items	Experimental conditions	Reaction accuracy rate	Spectral parameters
Staphylococcus aureus	66°C、 34min	30/32=93.8%	2128
Salmonella	66°C、 34min	30/32=93.8%	2135
Listeria monocytogenes	66°C、 34min	31/32=96.9%	2128
Hemorrhagic Escherichia coli	66°C、 34min	30/32=93.8%	2140
Escherichia coli (E. coli)	66°C、 34min	31/32=96.9%	2117

1.2 Feature Adjustments

1.3 Bug Fixes

2. Release note 002

Date: 2024.06.14

Changes from version 001

2.1 *New feature*

2.2 *Feature Adjustments*

2.3 *Bug Fixes*

- Fixed the risk of the peristaltic pump power potentially being set to zero after restoring factory settings.

File records

Version	Modification	Release
001	Initial release	2024.06.04
002	Fixed the risk of the peristaltic pump power potentially being zero.	2024.06.14